



MI08 Programmable Visible Spectrophotometer

NEW

www.camspec.co.uk

- ◆ 4 line x 20 character LCD display
- ◆ 9 button keypad
- ◆ Stepper motors for wavelength and filters
- ◆ Stores up to 200 methods/results
- ◆ Concentration calibration:
 - up to 8 standards
 - up to 3 replicates
 - linear or linear through zero
- ◆ RS232 serial output for printer or PC control
- ◆ Data output includes method number and name, wavelength, date and time, %T, Abs, concentration, unit
- ◆ Selectable wavelength calibration and dark current measurement
- ◆ 325 - 1000nm wavelength range
- ◆ 5nm bandwidth
- ◆ Low stray light (0.1%T at 340nm)
- ◆ Large sample compartment (up to 50mm pathlength cells)
- ◆ Large selection of accessories (details overleaf)
- ◆ Manual 4-cell changer, with a pair of 10mm glass cells included
- ◆ PC Application Software option for concentration curves, timescan and wavelength scanning

Accessory details overleaf >>>