



Spectronic News

SPECTRONIC ANALYTICAL INSTRUMENTS AN ISO 9001:2000 REGISTERED COMPANY

This newsletter gives you details on all of your accessory requirements for your Electrochemistry Meters.
Spectronic also offer bi-annual and annual maintenance contracts and calibration services.

ELECTROCHEMISTRY ACCESSORIES

SPECIAL OFFER

 * With effect from the 1st September 2008 until
 * 31st December 2008, Spectronic Analytical
 * Instruments are offering their full range of
 * electrodes and probes with a 10% discount along
 * with free postage and packing.
 *
 * To find out more about this offer visit our website
 * on www.spectronic.co.uk or contact our sales
 * team on 0113 2864536.

Electrodes & Probes

Spectronic Analytical Instruments are able to offer the following types of electrodes:
pH, Titrator, Conductivity, Dissolved oxygen, Ion selective, Redox and Reference electrodes.



For many of the leading manufacturers including:
ABB, Corning, Crison, Consort, Denver, Thermo Fisher, Hanna, Jenway, Merck, Metrohm, Mettler, Orion, Radiometer, Sartorius and WTW.
Combination Electrodes from £45
Food Industry Specific Electrodes also available

To discuss your requirements further contact our sales department on 0113 2864536 or to order on-line visit us @ Spectronic.co.uk

shop @ www.spectronic.co.uk

Ion Selective Electrodes



Available in the following:
Ammonium, Barium, Bromide, Cadmium, Calcium, Chloride, Cupric, Cyanide, Fluoride, Iodide, Lead, Mercury, Nitrate, Perchlorate, Potassium, Silver, Sodium, Sulphide, Thiocyanite and Water Hardness

- Main benefits and features include:
- No reference electrode needed.
 - No filling solution required.
 - Fully submersible and waterproof.
 - Solid state sensors.
 - Colour coded caps for identification in the field.
 - Can be left dry for long periods
 - Ideal for unskilled operatives.

Conductivity Probes



Available in the following:

- K=1 epoxy bodied
- K=1 glass
- K=1 glass F/T
- K=1 epoxy mini din
- K=0.1 epoxy bodied
- K=0.1 glass
- K=0.1 epoxy mini din
- K=10 epoxy bodied
- K=10 epoxy mini din

SPECTRONIC ANALYTICAL INSTRUMENTS
 TUDOR HOUSE, BARLEY HILL ROAD, GARFORTH, LEEDS, LS25 1DX
 Tel: 0113 286 4536 Fax: 0113 232 0424
 E-mail: sales@spectronic.co.uk www.spectronic.co.uk

9.2008
ELEC

Buffer Solutions & Capsules



	<u>Buffer Solutions</u>		<u>Buffer Capsules</u>
pH 4	500mls/5L	pH4	Box of 50
pH 7	500mls/5L	pH7	Box of 50
pH 9	500mls/5L	pH9	Box of 50
pH10	500mls/5L	pH10	Box of 50

Reference Filling Solutions

Electrode filling solution (KCl3.3m/AgNO3) (250ml)

KCl filling solution for D02 (250ml)

Electrode Storage & Cleaning Solution

Electrode cleaning solution (250ml)

Electrode storage solution (250ml)

Electrode storage bottle



Conductivity Standards

1413 μ S Conductivity standard (500ml)

12.88mS Conductivity standard (500ml)

111.9mS Conductivity standard (500ml)

10 μ S Conductivity standard (500ml)

84 μ S Conductivity standard (500ml)

TDS Standard 1382ppm

shop @ www.spectronic.co.uk

Spectronic Analytical Instruments have recently opened their on-line shop. Users can easily identify their consumable requirements and place a secure order on-line.



10% discount on orders over £100

To find out more about any of the accessories featured on this newsletter for your Electrochemistry Meter or any instrument in our range, log on to www.spectronic.co.uk

Spectronic accept payment via Visa, MasterCard, Visa Debit cards, Maestro, JCB cards and Visa company purchasing cards.